



OneTouch Verio[®] Pro+ blood glucose monitoring system

Record Book

ONETOUCH[®]
every touch is a step forward™

If this logbook is found,
please return to:

My Name _____

Surgery Name _____

CHP/LHB Name _____

Meter Serial Number _____

Customer Care Information

OneTouch® Customer Care



 **Call OneTouch Customer Care on
0800 121 200 (UK) or
1800 535 676 (Ireland).
Monday to Friday
9.00AM – 5.00PM**

 **CustomerCare@myOneTouch.co.uk**

 **www.myOneTouch.co.uk**

Before using the OneTouch Verio® Pro+ meter to do a control solution test, carefully read the Operator's Guide and inserts that come with the components of the system.

When testing multiple patients, it is important to follow your institution's policies and local standards and guidelines that apply to your region.


Control Solution For Your OneTouch Verio® Pro+ meter

Intended Use:

OneTouch Verio® Control Solution is used to check that the meter and test strips are working together properly and the test is performed correctly.

Please remember to complete all aspects of this diary as a record of your control testing.

What you need to perform a control solution test



- OneTouch Verio® Pro+ meter
- OneTouch Verio® test strip
- OneTouch Verio® Mid Control Solution

Performing a control solution test

- Whenever you open a new vial of test strips.
- Whenever using a new batch or lot of test strips.
- If you suspect the meter or test strips are not working properly (Contact Customer Care)
- If the meter has been dropped or damaged.
- If the meter has had repeated unexpected blood glucose results

Results outside the expected range may indicate:


- Failure to follow all the instructions for performing a control solution test described in your Meter Operator's Guide.
- Control solution is contaminated, expired or past its discard date
- Test strip or test strip vial is damaged, expired or past its discard date
- Meter, test strips and/or control solution were not all at the same temperature when the control solution test was performed
- A problem with the meter

If you continue to get control solution test results that fall outside the range printed on the control solution vial, Do Not use the meter, test strips, or control solution.

Contact OneTouch Customer Care on 0800 121 200 (UK) or 1800 535 676 (Ireland).

The control solution range printed on the control solution vial is for control solution tests only and is not a recommended range for your blood glucose level.

Should you experience any of the listed problems please call OneTouch Customer Care:

 **Call OneTouch Customer Care on
0800 121 200 (UK) or
1800 535 676 (Ireland).
Monday to Friday
9.00AM – 5.00PM**

Date	Time	Operator signature	Control Solution			Test Strip		Comments
			Expiry Date	Reading	Range	Expiry Date	Lot Number	

Notes:

Date	Time	Operator signature	Control Solution			Test Strip		Comments
			Expiry Date	Reading	Range	Expiry Date	Lot Number	

Notes:

Date	Time	Operator signature	Control Solution			Test Strip		Comments
			Expiry Date	Reading	Range	Expiry Date	Lot Number	

Notes:
