

OneTouch Reveal® Mobile Application Compatible Devices and System Requirements

To use the OneTouch Reveal® mobile application (app) you need:

- A compatible OneTouch® meter with Bluetooth® connectivity
- A compatible mobile device running a supported operating system
- To download the app to your mobile device from the App Store(SM) or Google Play

Compatible OneTouch® meter:

- OneTouch Verio Flex™

Compatible Mobile Devices & Operating Systems

The One Touch Reveal® mobile app is compatible with:

- Apple iPhone®, iPad® or iPod touch®, running iOS7 or later
- Google Android™ devices with BTLE (Bluetooth Low Energy) support and running Android Version 4.3, 4.4 and 5.1.

In order to use the Data sharing (Share) feature within the OneTouch Reveal® mobile app, your mobile device will need wifi or cellular connectivity to send by SMS message or email. See your mobile provider tariff for charges on SMS or data usage.

The Bluetooth® word mark and logos are registered trademarks owned by Bluetooth SIG, Inc., and any use of such marks by LifeScan Europe is under licence.

Apple, the Apple logo, iPhone, iPad, and iPod touch are trademarks of Apple Inc., registered in the U.S. and other countries. App Store is a service mark of Apple Inc.

iOS is a trademark of Apple Inc., registered in the U.S. and other countries.

Google, Google Play, Android and the Google Logo are registered trademarks of Google Inc.